



Liveo™ Pharma TPE Tubing

Pharmaceutical-grade thermoplastic elastomer (TPE) tubing dedicated to fluid transport and single-use bioprocessing applications.

APPLICATIONS

DuPont™ Liveo™ Pharma TPE Tubing is a sterilizable, weldable TPE tubing for biopharmaceutical processing for fluid transport and single-use bioprocessing applications:

- Fluid, media and solvent transport in biopharmaceutical processes and biotechnology
- Single-use assemblies
- Peristaltic pump applications
- Aseptic connection and disconnection without connectors (e.g., sampling)
- Applications where compatibility with industry benchmark TPE tubing is needed

FEATURES & BENEFITS

- Improved heat-welding, both to itself and to competitive offerings
- High tensile strength and burst resistance before and after welding
- Low extractables
- Good chemical resistance
- Minimal spallation after 24 hours of pumping
- Good clarity; stable clarity after sterilization
- Manufactured in ISO Class 7 cleanroom
- Comprehensive data package to accelerate qualification and validation under the form of a qualification guide



Liveo™ Pharma TPE Tubing

Common tubing product sizes and burst properties

Inner diameter (ID)		Outer diameter (OD)		Wall		Burst pressure			
						Autoclave ¹		Gamma ²	
inches	mm	inches	mm	inches	mm	psi	bar	psi	bar
1/8	3.2	1/4	6.4	1/16	1.6	94	6.5	89	6.1
3/8	9.5	5/8	15.9	1/8	3.2	87	6.0	85	5.9
1/2	12.7	3/4	19.1	1/8	3.2	71	4.9	70	4.8

Note: Burst pressure values are typical values not intended for writing specifications. Test method is based on ASTM D380-14.

¹Autoclave sterilization at 121°C for 30 minutes.

²Gamma treated 50 kGy dose.

SPECIAL NEEDS

DuPont has the capabilities to customize products to your specifications. Customization includes sizes, tolerances, cut lengths, bulk packaging and spooling. For easy identification, DuPont™ Liveo™ Pharma TPE Tubing also is available marked with product name and size. Please contact your DuPont representative to discuss your specific requirements.

Manufacturing environment

Liveo™ Pharma TPE Tubing is manufactured in the USA in an ISO Class 7 cleanroom at an ISO 13485:2016 certified site.



Liveo™ Pharma TPE Tubing

Characteristics

Delivery form	Tubing
Special characteristics	Translucent
Regulatory data package	USP <665> (Extractables and Leachables)* USP <85> (Endotoxins) USP <788> (Particulates) ISO 11737-1 (Bioburdens) ISO 10993 Biocompatibility tests (5: Cytotoxicity, 6: Muscle Implantation, 11: Acute Systemic Toxicity, 23: Intracutaneous Irritation) USP Class VI (Intracutaneous Reactivity, Acute Systemic Injection, Intramuscular Implant) USP <151> (Pyrogenicity) Elemental Impurities

*On autoclaved and gamma-radiated tubing.

Additional information

How to use	Information from “Chemical Resistance Guide for Elastomer IV,” published by Compass Publications, La Jolla, California, can be used as reference. It is, however, recommended that DuPont products be evaluated in the real conditions of use.
------------	--



Liveo™ Pharma TPE Tubing

Mechanical properties

Shore A hardness durometer	65 Shore A	ASTM D2240
----------------------------	------------	------------

Specific gravity	0.900	ASTM D792
------------------	-------	-----------

DuPont™ Liveo™ Pharma TPE Tubing ID 1/8" (3.2 mm) x OD 1/4" (6.4 mm)*

Stress at 200% strain, Die C	3.3 MPa (479 psi) ¹	ASTM D412
------------------------------	--------------------------------	-----------

	2.7 MPa (392 psi) ²	ASTM D412
--	--------------------------------	-----------

Maximum tensile strength, Die C	9.8 MPa (1,421 psi) ¹	ASTM D412
---------------------------------	----------------------------------	-----------

	8.2 MPa (1,189 psi) ²	ASTM D412
--	----------------------------------	-----------

Maximum elongation, Die C	902% ¹	ASTM D412
---------------------------	-------------------	-----------

	982% ²	ASTM D412
--	-------------------	-----------

DuPont™ Liveo™ Pharma TPE Tubing ID 3/8" (9.5 mm) x OD 5/8" (15.9 mm)

Stress at 200% strain, Die C	2.7 MPa (392 psi) ¹	ASTM D412
------------------------------	--------------------------------	-----------

	2.6 MPa (377 psi) ²	ASTM D412
--	--------------------------------	-----------

Maximum tensile strength, Die C	9.4 MPa (1,363 psi) ¹	ASTM D412
---------------------------------	----------------------------------	-----------

	7.4 MPa (1,073 psi) ²	ASTM D412
--	----------------------------------	-----------

Maximum elongation, Die C	1,134% ¹	ASTM D412
---------------------------	---------------------	-----------

	1,031% ²	ASTM D412
--	---------------------	-----------

*Tubing pulled apart as it is.

¹Autoclave sterilization at 121°C for 30 minutes.

²Gamma treated 50 kGy dose.



Liveo™ Pharma TPE Tubing

Mechanical properties (continued)

DuPont™ Liveo™ Pharma TPE Tubing ID 1/2" (12.7 mm) x OD 3/4" (19.1 mm)

Stress at 200% strain, Die C	2.7 MPa (392 psi) ¹	ASTM D412
	2.6 MPa (377 psi) ²	ASTM D412
Maximum tensile strength, Die C	9.0 MPa (1,305 psi) ¹	ASTM D412
	7.9 MPa (1,146 psi) ²	ASTM D412
Maximum elongation, Die C	1,128% ¹	ASTM D412
	1,073% ²	ASTM D412

¹Autoclave sterilization at 121°C for 30 minutes.

²Gamma treated 50 kGy dose.

Other properties

Shelf life 36 months¹

¹From the date of manufacture when stored in the original unopened containers.

Sterilization methods

- Autoclave steam, 30 minutes at 121°C
- Gamma radiation, 50 kGy

Revised: 2021-10-07

Page: 5 of 5

dupont.com

The information set forth herein is furnished free of charge, is based on technical data that DuPont believes to be reliable, and represents typical values that fall within the normal range of properties. This information relates only to the specific material designated and may not be valid for such material used in combination with other materials or in other processes. It is intended for use by persons having technical skill, at their own discretion and risk. This information should not be used to establish specification limits nor used alone as the basis of design. Handling precaution information is given with the understanding that those using it will satisfy themselves that their particular conditions of use present no health or safety hazards and comply with applicable law. Since conditions of product use and disposal are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information. As with any product, evaluation under end-use conditions prior to specification is essential. Nothing herein is to be taken as a license to operate or a recommendation to infringe on patents.

CAUTION: Do not use DuPont materials in medical applications involving implantation in the human body or contact with internal body fluids or tissues unless the material has been provided from DuPont under a written contract or other acknowledgement that is consistent with the DuPont policy regarding medical applications and expressly acknowledges the contemplated use. For further information, please contact your DuPont representative.

DuPont's sole warranty is that our products will meet our standard sales specifications in effect at the time of shipment. Your exclusive remedy for breach of such warranty is limited to refund of purchase price or replacement of any product shown to be other than as warranted. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, DUPONT SPECIFICALLY DISCLAIMS ANY OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY, OR NON-INFRINGEMENT. DUPONT DISCLAIMS LIABILITY FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES.

DuPont™, the DuPont Oval Logo, and all trademarks and service marks denoted with ™, SM or ® are owned by affiliates of DuPont de Nemours, Inc. unless otherwise noted.
© 2021 DuPont. All rights reserved.

AGP16478M-22