OUPONT>

Liveo™ Dimethiconol Blend 20

Silicone Topical Excipients

Excipient for pharmaceutical applications

Liveo[™] Dimethiconol Blend 20 is a blend of an ultra high molecular weight hydroxyterminated polydimethylsiloxane (dimethiconol) in a low viscosity and non-volatile silicone fluid (dimethicone, 20 cSt).

Composition

INCI name: Dimethicone and Dimethiconol, CAS number: 70131-67-8 and 63148-62-9

Applications, Features & Benefits:

Liveo[™] Dimethiconol Blend 20 is commonly used in pharmaceutical formulations thanks to its long lasting, film barrier and wash-off resistance properties. Such silicone gum blends are currently used in a range of dermatological and pharmaceutical applications such as protective and barrier creams, gels for scar treatment and formulations for skin and scalp treatments. Clear, Non-volatile fluid as carrier, Substantive, Non-occlusive, Film former, Long-lasting, Wash off resistance, Smooth and lubricious feel

Packaging:

Liveo[™] Dimethiconol Blend 20 is available in pails and in drums. Samples for Liveo[™] Dimethiconol Blend 20 are available in bottles.

Product information

Colour Polydimethylsiloxane content	Clear - 12 %	PhEur 1470
Rheological properties Viscosity	387 mPa.s	
Thermal properties Flash point	200 °C	ASTM D 92
Other properties Density Refractive index	952 kg/m³ 1.401 -	ISO 1183
Storage and stability Shelf life	24 months	
Additional Information		
How to use	Liveo™ Dimethiconol Blend 20 is an easy to use blend of dimethicone and dimethiconol. It is compatible with a range of topical excipients. Liveo™ Dimethiconol Blend 20 can be added to various anhydrous or water	

based formulations. Liveo[™] Dimethiconol Blend 20 can be used alone or blended with other excipient to provide a fluid base for a variety of



Liveo™ Dimethiconol Blend 20

Silicone Topical Excipients

formulation and provide a smooth, long lasting film on skin. The presence of Liveo[™] Dimethiconol Blend 20 can increase the viscosity of water-in-oil emulsions.

When preparing an emulsion containing Liveo[™] Dimethiconol Blend 20, it is recommended to add the oil(s) to the silicone blend under high shear speed; the emulsion can then potentially be processed through a homogenizer to improve its stability. The blend is easily emulsified in water with standard emulsifiers and normal emulsification techniques.

The information set forth herein is furnished free of charge, is based on technical data that DuPont believes to be reliable, and represents typical values that fall within the normal range of properties. This information relates only to the specific material designated and may not be valid for such material used in combination with other materials or in other processes. It is intended for use by persons having technical skill, at their own discretion and risk. This information hould not be used to establish specification limits nor used alone as the basis of design. Handling precaution information is given with the understanding that those using it will satisfy themselves that their particular conditions of use present no health or safety hazards and comply with applicable law. Since conditions of product use and disposal are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information. As with any product, evaluation under end-use conditions prior to specification is essential. Nothing herein is to be taken as a license to operate or a recommendation to infringe on patents.

CAUTION: Do not use DuPont materials in medical applications involving implantation in the human body or contact with internal body fluids or tissues unless the material has been provided from DuPont under a written contract or other acknowledgement that is consistent with the DuPont policy regarding medical applications and expressly acknowledges the contemplated use. For further information, please contact your DuPont representative.

DuPont's sole warranty is that our products will meet our standard sales specifications in effect at the time of shipment. Your exclusive remedy for breach of such warranty is limited to refund of purchase price or replacement of any product shown to be other than as warranted. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, DUPONT SPECIFICALLY DISCLAIMS ANY OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY, OR NON-INFRINGEMENT. DUPONT DISCLAIMS LIABILITY FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES.

DuPont[™], the DuPont Oval Logo, and all trademarks and service marks denoted with [™], SM or [®] are owned by affiliates of DuPont de Nemours, Inc. unless otherwise noted. © 2021 DuPont. All rights reserved.